

ORP/PH tester/Simulator/Model 207

FUNCTIONS

For a calibration of PH & ORP meter & for checking of sensor response.

FEATURES

- Portable, ideal for laboratory & field calibration of pH & ORP systems.
- Universal Calibrator.
- Temp Calibration (PT 1000/NTC 300Ω/PT 100)
Range :0⁰-100⁰ C
- Multipoint calibration of pH, mV & temperature.
- Millivolt Calibration (Range: -1000mv to +1000mv)
- PH Calibration (0-14 PH)
- Display: (2.5" X 1")
- PH/ ORP Probe health check for routine maintenance.
- Rugged Aluminum body.



Ordering Information

207 - pH Tester

T1 - 300 Ω NTC
T2 - PT 1000
T3 - PT 100